




Polypropylene Cell Spreader LD-PCS101L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polypropylene Cell Spreader LD-PCS101L

Features:

- Material: Polypropylene
- Ideal for bacterial CFU counts
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

Specifications:

Size	21.5 x 3.0 mm
Color	Blue
Shape	Triangular
Packaging	1 pcs/bag
Sterilization	Gamma ray sterilized
Quantity per case	800