



ABS Centrifuge Tube Rack LD-CTR102L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

ABS Centrifuge Tube Rack LD-CTR102L

Features:

- Temperature variance range (-80 Â°C to 121 Â°C)
- Material: Acrylonitrile Butadiene Styrene
- Autoclavable
- Availability of color customization

Specifications:

Packaging	Plastic Bag
No. of wells	20 + 30
Tube Capacity	15 mL & 50 mL (Multi use)
Quantity per case	50