



Conical Bottom Centrifuge Tubes LD-CCT201L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Conical Bottom Centrifuge Tubes LD-CCT201L

Features:

- Maximum RCF 12000 x g
- Material: Virgin Polypropylene (PP)
- Working Temperature is -80 °C to 121 °C
- Transparent wall permits easy observation of tube contents
- Graduation in white ink
- Frosted marking area
- Flat-top screw cap for free handling

Specifications:

Capacity	50 mL
Packaging	25 pcs/Foam Rack
Length of tube	11.7 cm
Diameter of Cap	34 mm
Quantity per Case	500
Diameter outside tube	30 mm