



Microcentrifuge Tubes LD-MCT201L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Microcentrifuge Tubes LD-MCT201L

Features:

- Maximum RCF 28000 x g
- Working Temperature is -90 °C to 121 °C
- Conical bottom design
- Autoclavable temperature is +121 °C
- Material: Polypropylene
- Frosted marking area on tube and cap
- Attached cap

Specifications:

Type	Ultra Low Retention
Capacity	1.5/1.6/1.7 mL
Packaging	500 pcs/Bag
Length of tube	41 mm
Quantity per Case	5000
Diameter of the tube	11 mm