



Plastic Microcentrifuge Tube Rack LD-MTR100L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Plastic Microcentrifuge Tube Rack LD-MTR100L

Features:

- Easy removal and insertion of tubes
- It can be autoclaved
- It is compatible with 1.5 mL, 2.0 mL, 5.0 mL and 15.0 mL
- Available in two colors (Orange and Blue) for better color coding

Specifications:

Places	50 wells (10 x 2 array)
Capacity	1.5/ 2.0 mL
Packaging	Individual
Diameter of well	12.2 mm
Dimension (L*W)	203*122 mm
Quantity per Case	50