



ABS Cryovials Holder LD-CVH100L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

ABS Cryovials Holder LD-CVH100L

Features:

- Specially designed for cryovials
- Autoclavable
- Chemically resistant to alcohols and mild organic solvents
- Temperature range is stable from -196 °C to 121 °C

Specifications:

Capacity	2.0 mL
Diameter	13 mm
Material	Acrylonitrile butadiene styrene
Packaging	10 pcs/Bag
Number of Wells	25
Quantity per Case	100
Dimension (L*W*H)	76*76*54 mm