




ABS Cryovials Holder LD-CVH102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

ABS Cryovials Holder LD-CVH102L

Features:

- Specially designed for cryovials
- Autoclavable
- Chemically resistant to alcohols and mild organic solvents
- Temperature range is stable from -196 °C to 121 °C

Specifications:

Capacity	2.0 mL
Diameter	11 mm
Material	Acrylonitrile butadiene styrene
Packaging	10 pcs/Bag
Number of Wells	100
Quantity per Case	50
Dimension (L*W*H)	132*132*54 mm