



External Threaded Cryo Vials LD-ECV100L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

External Threaded Cryo Vials LD-ECV100L

Features:

- Equipped with Medical class silicon seal-ring
- Working temperature is -196 °C to 121 °C
- DNase, RNase and Pyrogenic free
- Pre-Sterilized
- Self-standing but Conical bottom
- Used for storing biological material, human or animal cells

Specifications:

Type	Self-standing
Capacity	1.0 mL
Packaging	25 pcs/Bag
Length of Tube	35.50 mm
Diameter of Tube	12.50 mm
Quantity per Case	500