



## **External Threaded Cryo Vials LD-ECV104L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# External Threaded Cryo Vials LD-ECV104L

## Features:

- Equipped with Medical class silicon seal-ring
- Working temperature is -196 °C to 121 °C
- DNase, RNase and Pyrogenic free
- Pre-Sterilized
- Self-standing but Conical bottom
- Used for storing biological material, human or animal cells

## Specifications:

|                   |               |
|-------------------|---------------|
| Type              | Self-standing |
| Capacity          | 5.0 mL        |
| Packaging         | 25 pcs/Bag    |
| Length of Tube    | 90.50 mm      |
| Diameter of Tube  | 12.50 mm      |
| Quantity per Case | 500           |