



Internal Threaded Cryo Vials LD-ICV101L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Internal Threaded Cryo Vials LD-ICV101L

Features:

- Working temperature is -196 °C to 121 °C
- Caps are made of Medical class Silicon
- DNase, RNase and Pyrogenic free
- Pre-Sterilized
- Self-standing but Conical bottom
- Used for storing biological material, human and animal cells

Specifications:

Type	Self-standing
Capacity	1.8 mL
Packaging	25 pcs/Bag
Length of Tube	48.3 mm
Diameter of Tube	12.40 mm
Quantity per Case	500