



PCR Well Plates LD-PWP200L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

PCR Well Plates LD-PWP200L

Features:

- Maximum G-Force is 3500 RCF
- Material: Highest Grade Polypropylene
- Compatible with range of PCR and qPCR machines
- Thin-walled plates for optimum PCR and real-time (RT, qPCR) results
- Pre-sterilized
- RNase, DNase and pyrogenic free

Specifications:

Type	Semi-Skirt
Color	Nature
Capacity	0.1 mL
Packaging	15 pcs/box
Quantity per Case	150