



PCR Tube LD-PTB104L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

PCR Tube LD-PTB104L

Features:

- Material: Polypropylene
- Ultra-thin walls for efficient thermal transfer and maximum yield
- Domed or flat insertion caps ensure a tight fit that prevents sample evaporation
- Optically clear cap for real time qPCR

Specifications:

Cap Type	Flat Cap
Packaging	200 pcs/Bag, 20 Bags/Case
Specification	0.2 mL, 12 Strip
Quantity per case	4000