



Racked Bar-Coded Cryogenic Vials LD-BCV106L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Racked Bar-Coded Cryogenic Vials LD-BCV106L

Features:

- Material: Polypropylene
- Low Temperature Storage: -196 Â°C to 121 Â°C
- Sterilization Level: SAL 10-6
- Chemically resistant and scratch-proof 2D Laser Etched code on Bottom
- U-shaped for maximum working volume
- Barcoded rack for tracking
- Inner Diameter: 0.7 mm
- Cap Color: Transparent, White, Blue, Orange and Green
- DNase, RNase and Pyrogenic free

Specifications:

Height	33.3 mm
Barcode	1D
Cap Type	Screw Cap
Capacity	1.4 mL
Packaging	96 pcs/Rack, 48 Racks/Case
Working Volume	1.15 mL
Quantity per Case	4608