



Fixed Volume Fully Autoclavable Pipettes LD-FVP110L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Fixed Volume Fully Autoclavable Pipettes LD-FVP110L

Features:

- Autoclavable at 121 Â°C
- Material: Polycarbonate
- Excellent durability and chemical resistance
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

Specifications:

Test Volume	2000 µL
Accuracy error	0.3 %
Precision error	0.15 %