



Fixed Volume Pipettes LD-FVP200L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Fixed Volume Pipettes LD-FVP200L

Features:

- Autoclavable at 121 Â°C
- Material: Polycarbonate
- Excellent durability and chemical resistance
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

Specifications:

Test Volume	5 µL
Accuracy error	1.3 %
Precision error	1.2 %