



Fixed Volume Pipettes LD-FVP207L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Fixed Volume Pipettes LD-FVP207L

Features:

- Autoclavable at 121 Â°C
- Material: Polycarbonate
- Excellent durability and chemical resistance
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

Specifications:

Test Volume	250 ?L
Accuracy error	0.4 %
Precision error	0.2 %