



Variable Volume Fully Autoclavable Pipettes LD-VVP106L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Variable Volume Fully Autoclavable Pipettes LD-VVP106L

Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Available in 0.1 µL to 5000 µL volume
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

Specifications:

Increment	1 µL
Test Volume	200 µL <hr> 100 µL <hr> 50 µL
Volume Range	50 - 200µL
Accuracy error	0.60 % <hr> 0.80 % <hr> 1.00 %
Precision error	0.15 % <hr> 0.30 % <hr> 0.40 %