



Variable Volume Pipette LD-VVP208L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Variable Volume Pipette LD-VVP208L

Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Available in 0.1 µL to 10000 µL volume
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

Specifications:

Increment	5 µL
Test Volume	1000 µL <hr> 500 µL <hr> 200 µL
Volume Range	200 - 1000µL
Accuracy error	0.60 % <hr> 0.70 % <hr> 0.90 %
Precision error	0.20 % <hr> 0.25 % <hr> 0.30 %