



Pipette Controller LD-PPC101L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Pipette Controller LD-PPC101L

Features:

- Compatible with Glass and Plastic Pipette from 0.1 -100mL
- Single Handed Operation
- Eight modes speed selection for aspiration and dispensing different liquids
- LCD display for showing low battery warning and speed settings
- Battery Service Life
- High capacity Li-ion battery enabling long runtime of operation
- 0.45µm replaceable hydrophobic filter
- Rechargeable during use
- Battery life lasts for more than 8 hours of Intermittence use

Specifications:

Filter	0.45 µm Hydrophobic
Weight	200g
Charging Time	2 to 3 Hours
Aspiration Speed	25 mL < 5 s (8 shift)
Dispensing Speed	Motor (8 shift) / Gravity
Quantity Per case	1