



Filter Pipette Tips LD-FPT115L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Filter Pipette Tips LD-FPT115L

Features:

- Universal compatibility
- Material: Polypropylene
- Micro point design
- Pre sterilized
- RNase, DNase and Pyrogenic free
- Optional Color “ Yellow (FPT112L and FPT113L)
- Optional Color “ Blue (FPT114L, FPT115L, FPT116L and FPT117L)

Specifications:

Type	Long
Length	86.4 mm
Capacity	1000 µL
Packaging	1000 Tips/Bag, 10 Bags/Case
Quantity per Case	10000