



Ultra Low retention Pipette Tips LD- LRT117L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Ultra Low retention Pipette Tips LD-LRT117L

Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121Â°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color “ Yellow (LRT106L, LRT107L, LRT108L, LRT109L, LRT110L and LRT111L)
- Optional Color “ Blue (LRT112L, LRT113L, LRT114L, LRT115L, LRT116L and LRT117L)

Specifications:

Type	Long
Capacity	1250/1500 µL
Packaging	768 Tips/Refill unit , 10 Units/Case
Length of tube	103.20 mm
Quantity per Case	7680