



Pipette Tip for Eppendorf LD-SPT102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Pipette Tip for Eppendorf LD-SPT102L

Features:

- Material: Polypropylene
- Autoclavable
- Smooth inner wall to avoid leakage and sample residual
- Pre Sterilized
- DNase, RNase and pyrogen free

Specifications:

Type	Non Filter
Color	Yellow
Capacity	200 ?L
Packaging	1000 pcs/Pack, 30 Packs/Case
Compatibility	Eppendorf
Quantity Per Case	30000