



Pipette Tip for Finn LD-SPT108L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Pipette Tip for Finn LD-SPT108L

Features:

- Material: Polypropylene
- Autoclavable
- Smooth inner wall to avoid leakage and sample residual
- Pre Sterilized
- DNase, RNase and pyrogen free

Specifications:

Type	Non-Filter
Color	Transparent
Capacity	200 µL
Packaging	1000 pcs/Pack, 30 Packs/Case
Compatibility	Finn
Quantity Per Case	30000