



Pipette Tip for Finn LD-SPT112L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Pipette Tip for Finn LD-SPT112L

Features:

- Material: Polypropylene
- Autoclavable
- Smooth inner wall to avoid leakage and sample residual
- Pre Sterilized
- DNase, RNase and pyrogen free

Specifications:

Type	Non-Filter
Color	Transparent
Capacity	5 mL
Packaging	300 pcs/Pack, 10 Packs/Case
Compatibility	Finn
Quantity Per Case	3000