



## **Pipette Tip for Gilson LD-SPT114L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Pipette Tip for Gilson LD-SPT114L

## Features:

- Material: Polypropylene
- Autoclavable
- Smooth inner wall to avoid leakage and sample residual
- Pre Sterilized
- DNase, RNase and pyrogen free

## Specifications:

|                   |                              |
|-------------------|------------------------------|
| Type              | Non-Filter                   |
| Color             | Yellow                       |
| Capacity          | 200 ?L                       |
| Packaging         | 1000 pcs/Pack, 30 Packs/Case |
| Compatibility     | Gilson                       |
| Quantity Per Case | 30000                        |