




Universal Pipette Tips LD-UPT105L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Universal Pipette Tips LD-UPT105L

Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121Â°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color “ Yellow (UPT106L, UPT107L, UPT108L, UPT109L, UPT110L and UPT111L)
- Optional Color “ Blue (UPT112L, UPT113L, UPT114L, UPT115L, UPT116L, UPT117L and UPT118L)

Specifications:

Type	Short
Capacity	10 µL
Packaging	960 Tips/Refill Unit, 10 Unit/Case
Length of tube	31.70 mm
Quantity per Case	9600