




Universal Pipette Tips LD-UPT107L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Universal Pipette Tips LD-UPT107L

Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121Â°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color “ Yellow (UPT106L, UPT107L, UPT108L, UPT109L, UPT110L and UPT111L)
- Optional Color “ Blue (UPT112L, UPT113L, UPT114L, UPT115L, UPT116L, UPT117L and UPT118L)

Specifications:

Capacity	200 µL
Packaging	1000 Tips/Bag, 10 Bags/Case
Length of tube	50 mm
Quantity per Case	10000