



Polystyrene Serological Pipette LD- PSP104L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polystyrene Serological Pipette LD-PSP104L

Features:

- Good optical Transmittance
- Material: Polystyrene
- Gamma ray sterilized
- Equipped with Porex filter
- DNase, RNase and endotoxin free

Specifications:

Capacity	25 mL
Packaging	200 Pcs/Case
Color of Mark	Purple
Min. graduation	2/10 mL
Quantity per Case	200