



PolyEthylene Transfer Pipette LD-PTP100L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

PolyEthylene Transfer Pipette LD-PTP100L

Features:

- Material: Low-density Polyethylene
- Heat seal on the tip can prevent liquid leakage
- Gamma ray sterilized
- Non-toxic
- High Transparency and Elasticity

Specifications:

Length	143 mm
Capacity	1 ml
Drop Vol.	42 ?L
Packaging	500 pcs/pk
Quantity Per Case	5000