



PolyEthylene Transfer Pipette LD-PTP101L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

PolyEthylene Transfer Pipette LD-PTP101L

Features:

- Material: Low-density Polyethylene
- Heat seal on the tip can prevent liquid leakage
- Gamma ray sterilized
- Non-toxic
- High Transparency and Elasticity

Specifications:

Length	145 mm
Capacity	1 mL
Drop Vol.	45 ?L
Packaging	500 pcs/pk
Quantity Per Case	5000