



Plastic Pipette Stand LD-PPS101L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Plastic Pipette Stand LD-PPS101L

Features:

- Material: Polymethyl-methacrylate and Acrylonitrile butadiene styrene
- Available in 3 different models
- Ensures easy removal and replacement of pipettes with weighted base
- Convenient and safe storage

Specifications:

Type	Z Type
Color	White
Holds	6
Material	PMMA
Packaging	Bulk Packing
Quantity Per Case	10