



Plastic Pipette Stand LD-PPS201L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Plastic Pipette Stand LD-PPS201L

Features:

- Material: Polymethyl-methacrylate and Acrylonitrile butadiene styrene
- Available in 3 different models
- Ensures easy removal and replacement of pipettes with weighted base
- Convenient and safe storage

Specifications:

Type	Disc Type
Color	White
Holds	6
Material	ABS
Packaging	1 pc/Box, 20 Box/Case
Quantity Per Case	20